

Protocol for Freezing Cells

- 1) Grow cells until full but not completely confluent
- 2) Harvest cells into 15mL conical tube
- 3) Spin cells at 1,000 RPM (173xg) for 5min.
- 4) Aspirate media and resuspend in 1mL freezing media per freezing vial
- 5) Transfer vials to -80C (use a Styrofoam rack if available)
- 6) For long-term storage transfer tubes from -80C into liquid nitrogen

10cm plate freeze into 1-2 vials

15cm plate freeze into 3-4 vials

Freezing Media

20mL FBS

25mL media

5mL DMSO

Filter sterilize using 50mL Steriflip vacuum filter unit

Thawing cells

- 1) Remove vial from -80C or liquid nitrogen and place in 37C water bath just long enough for thaw.
- 2) Immediately upon thaw, transfer thawed cells into a 15mL conical tube containing 5mL fresh cell culture media
- 3) Spin cells at 1,000 RPM (173xg) for 5min.
- 4) Aspirate cell culture media and resuspend in fresh cell culture media
- 5) Plate cells
- 6) Next day: remove media from plate and add fresh media or split cells if the plate is full